

DNA NEW TABLE SETTINGS IN ID-X

A.1 **TOOLS>GENEMAPPER ID-X MANAGER (CTRL+M)**

A.2 Select the **TABLE SETTINGS TAB**

A.3 Select **NEW...** at the bottom of the dialog box

A.4 In the **GENERAL TAB:**

A.4.1 Name: name table setting

A.4.2 Security group: Genemapper *ID-X* security group

A.4.3 Description: leave blank

A.5 In the **SAMPLES TAB:**

A.5.1 Column Settings: select which features will be shown in the table

A.5.2 Font Settings: set desired font

A.5.3 Sort By: set to analysts preference

A.6 In the **GENOTYPES TAB:**

A.6.1 Column Settings: select which features will be shown in the table

A.6.2 Font Settings: set desired font

A.6.3 Sort By: set to analysts preference

A.6.4 Show number of alleles: 2

A.6.5 Allele position: unchecked

A.7 Click **OK** to **SAVE**

Document ID	Revision	Approval	Date Published
1813	8	Supervising Criminalist - Biology	12/28/2018 10:36:47 AM